



Email : sales@naugramedical.com

Phone: 9896600003

Product Name :
Semi Fowler Steel Bed Manual

Product Code :
NMHOPB00029



Description :

Semi Fowler Steel Bed Manual India China

Technical Specification :

Semi Fowler Steel Bed Manual with mechanically operated back rest by super-smooth crank mechanism. Size: L 2120mm

W 940mm

H 500mm. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Epoxy coated mild steel frame work and 2 section top. Mounted on protective stumps. Provision for I.V. Rod on both sides of the bed. Freight saving knock-down construction



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,