



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Semi Fowler Bed Manual

Product Code :
NMHOPB00030



Description :

Semi Fowler Bed Manual

Technical Specification :

Semi Fowler Bed Manual Bed with mechanically operated back rest by super-smooth crank mechanism. Size: L 2120mm

W 940mm

H 500mm. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy) Epoxy coated mild steel frame work and 2 section perforated top. Mounted on protective stumps. Provision for I.V. Rod on both sides of the bed. Freight saving knock-down construction



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,