



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Standard Semi-Fowler Bed

**Product Code :**  
NMHOPB00031



**Description :**

Standard Semi-Fowler Bed Standard semi-fowler bed with mechanically operated back rest by crank mechanism. Size: 2030  $\times$  900  $\times$  530 mm. Epoxy coated mild steel frame work and 2 section top. Tubular head & foot bows made of stainless steel. Mounted on protective stumps. Provision for I.V. Rod. Freight saving knock-down construction



**Medical Hospital Equipments Manufacturers**

**Website:** Medical Hospital Equipments Manufacturers,