



Email : sales@naugramedical.com

Phone: 9896600003

**Product Name :**  
Hospital Fowler Bed Manual

**Product Code :**  
NMHOPB00040



**Description :**

**Hospital Fowler Bed Manual India China**

**Technical Specification :**

Hospital Fowler Bed Manual with mechanically operated back rest and knee rest by super-smooth crank mechanism. Size: L 2140mm

W 940mm

H 500mm. Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Equipped with roller Bumper. Epoxy coated mild steel frame work and 4 section perforated top. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod



**Medical Hospital Equipments Manufacturers**

Website: [SITE](http://SITE)