



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Naugramedical steel Fowler Bed Manual

Product Code :
NMHOPB00042



Description :

Naugramedical steel Fowler Bed Manual India China

Technical Specification :

Naugramedical steel Fowler Bed Manual with mechanically operated back rest and knee rest by super-smooth crank mechanism. Size: L 2140mm W 940mm

H 530mm. Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Epoxy coated mild steel frame work and 4 section perforated top. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,