



Email : sales@naugramedical.com

Phone: 9896600003

Product Name :
Fowler Knee Rest Bed Manual

Product Code :
NMHOPB00043



Description :

Fowler Knee Rest Bed Manual India China

Technical Specification :

Fowler Knee Rest Bed Manual with mechanically operated back rest and knee rest by super-smooth crank mechanism.

Size: L 2190mm

W 940mm

H 500mm. Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Epoxy coated mild steel frame work and 4 section perforated top. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod



Medical Hospital Equipments Manufacturers

Website: SITE

