

Product Name :
Naugramedical Fowler Bed Manual 2 Function

Product Code :
NMHOPB00044



Description :

Naugramedical Fowler Bed Manual 2 Function India China

Technical Specification :

Naugramedical Fowler Bed Manual 2 Function with mechanically operated back rest and knee rest by super-smooth crank mechanism. Size: L 2140mm W 940mm

H 500mm. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy) Epoxy coated mild steel frame work and 4 section perforated top. 125mm dia noiseless castors (all with individual brakes) Provision for I.V. Rod