



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Hospital Fowler Bed Manual 2 Function

Product Code :
NMHOPB00046



Description :

Hospital Fowler Bed Manual 2 Function India China

Technical Specification :

Hospital Fowler Bed Manual 2 Function with mechanically operated back rest and knee rest by super-smooth crank mechanism. Size: L 2190mm

W 980mm

H 500mm. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy) Epoxy coated mild steel frame work and 4 section perforated top. 125mm dia noiseless castors (all with individual brakes) Provision for I.V. Rod on both sides of the bed. Freight saving knock-down construction



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,