



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Hospital Fowler Bed Manual 2 Function

**Product Code :**  
NMHOPB00046



**Description :**

**Hospital Fowler Bed Manual 2 Function India China**

**Technical Specification :**

Hospital Fowler Bed Manual 2 Function with mechanically operated back rest and knee rest by super-smooth crank mechanism. Size: L 2190mm

W 980mm

H 500mm. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy) Epoxy coated mild steel frame work and 4 section perforated top. 125mm dia noiseless castors (all with individual brakes) Provision for I.V. Rod on both sides of the bed. Freight saving knock-down construction



**Medical Hospital Equipments Manufacturers**

Website: [SITE](http://SITE)