



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Fowler Bed Manual 2 Function

**Product Code :**  
NMHOPB00048



**Description :**

Fowler Bed Manual 2 Function

**Technical Specification :**

Fowler Bed Manual 2 Function with mechanically operated back rest and knee rest by super-smooth crank mechanism. Size: L 2140mm

W 940mm

H 500mm. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. D type Design. Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy) Epoxy coated mild steel frame work and 4 section perforated top. 125mm dia noiseless castors (all with individual brakes) Provision for I.V. Rod on both sides of the bed. Freight saving knock-down construction



**Medical Hospital Equipments Manufacturers**

Website: [SITE](http://SITE)