



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Fowler Back Rest Bed Manual 2 Function

**Product Code :**  
NMHOPB00050



**Description :**

Fowler Back Rest Bed Manual 2 Function

**Technical Specification :**

Fowler Back Rest Bed Manual 2 Function with mechanically operated back rest and knee rest by super-smooth crank mechanism. Size: L 2140mm

W 980mm

H 500mm. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy) Epoxy coated mild steel frame work and 4 section perforated top. Mounted on protective stumps. Provision for I.V. Rod on both sides of the bed. Freight saving knock-down construction



**Medical Hospital Equipments Manufacturers**

Website: [SITE](http://SITE)