



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Standard Fowler Knee Rest Bed

Product Code :
NMHOPB00052



Description :

Standard Fowler Knee Rest Bed

Technical Specification :

Standard Fowler Bed Standard fowler bed with mechanically operated back rest and knee rest by crank mechanism.

Size: 2030 900 500 mm.

Epoxy coated mild steel frame work and 4 section top. Stainless steel head & foot panels. Four swivel castor (all with individual brakes). Provision for I.V. Rod. Freight saving knock-down construction



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,