



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Standard Fowler Bed

Product Code :
NMHOPB00053



Description :

Standard Fowler Bed Standard fowler bed with mechanically operated back rest and knee rest by crank mechanism. Size: 2030 \times 900 \times 530 mm. Epoxy coated mild steel frame work and 4 section top. Tubular head & foot bows made of stainless steel. Mounted on protective stumps. Provision for I.V. Rod. Freight saving knock-down construction



Medical Hospital Equipments Manufacturers

Website: SITE