



Email : sales@naugramedical.com

Phone: 9896600003

Product Name :
Standard Fowler Steel Bed

Product Code :
NMHOPB00054



Description :

Standard Fowler Steel Bed

Technical Specification :

Standard Fowler Bed Standard fowler bed with mechanically operated back rest and knee rest by crank mechanism.

Size: 2030 900 500 mm.

Epoxy coated mild steel frame work & 4 section top.. Mounted on protective stumps. Provision for I.V. Rod. Freight saving knock-down construction



Medical Hospital Equipments Manufacturers

Website: [SITE](#)