



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Standard Fowler Steel Tubes Bed

Product Code :
NMHOPB00055



Description :

Standard Fowler Steel Tubes Bed

Technical Specification :

Standard Fowler Bed Size 206L x 90W x 60H cms Frame work made of strong & precise steel tubes. Epoxy Powder Coated Finish. 4 Section Top made of perforated CRCA Sheet with mattress guards. Back Rest and Knee Rest Positions manoeuvred by Separate screw from Foot End. Tubular head and foot bows of unequal height with vertical tube supports. Provision for four saline locations. Bed Mounted on 100mm Dia Castor $1\frac{1}{2}$ two with brakes. Knockdown Construction



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,