



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Fowler Bed Manual Two Function

Product Code :
NMHOPB00059



Description :

Fowler Bed Manual Two Function

Technical Specification :

Fowler Bed Manual Two Function with mechanically operated back rest and knee rest by super-smooth crank mechanism.

Size: L 2120mm W 940mm

H 500mm. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Epoxy coated mild steel frame work and 4 section perforated top. Mounted on protective stumps. Provision for I.V. Rod on both sides of the bed. Freight saving knock-down construction



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,