



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Fowler Bed, Manual (with Provision for Bed Pan)

Product Code :
NMHOPB00062



Description :

Fowler Bed, Manual (with Provision for Bed Pan) Fowler Bed with mechanically operated back rest 0 to 75° and knee rest 0 to 40° by crank mechanism. Size: 2000 x 900 x 500 mm. Provision for Bed pan. With side railings. Epoxy coated mild steel frame work and 4 section top. Mounted on protective stumps. Provision for I.V. Rod on both sides of the bed. Freight saving knock-down construction



Medical Hospital Equipments Manufacturers

Website: SITE