

Product Name :

Infant Bed / Child Cot / Baby Bassinet (with Trendelenburg & Reverse Trendelenburg 12½°)

Product Code :

BCIVP10001

**Description :**

Infant Bed / Child Cot / Baby Bassinet (with Trendelenburg & Reverse Trendelenburg 12½°) Size: 870mm L x 530mm W x 920-1070mm H. Height Adjustable by Gas Spring 920-1070mm. Moulded Plastic Crib (without joints) with clear visibility for observer. Smooth Rotating Castors 100mm Dia will provide silent transportation. Brakes on opposite castors. With 12½° Trendelenburg (TR) & Reverse Trendelenburg (RTR) position