

**Product Name :**  
Infant Bed / Child Cot / Baby Bassinet (with Trendelenburg  
& Reverse Trendelenburg 12½°)

**Product Code :**  
BCIVP10001



**Description :**

Infant Bed / Child Cot / Baby Bassinet (with Trendelenburg & Reverse Trendelenburg 12½°) Size: 870mm L  
530mm W 920-1070mm H. Height Adjustable by Gas Spring 920-1070mm. Moulded Plastic Crib  
(without joints) with clear visibility for observer. Smooth Rotating Castors 100mm Dia will provide silent  
transportation. Brakes on opposite castors. With 12½° Trendelenburg (TR) & Reverse Trendelenburg (RTR)  
position