

Product Name :
Infant Bed / Child Cot with Plastic Moulded Crib

Product Code :
BCIVP10007



Description :

Infant Bed / Child Cot with Plastic Moulded Crib Size: 810 × 460 × 810 mm. Removable moulded Plastic Crib (without joints). With Trendelenburg (TR) Position. (Inclination upto 12°) Framework made of stainless steel. 1 S.S. shelf at the bottom. 4 swivel castor (2 with brakes)

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,