

Product Name :
Baby Bassinet Deluxe (With poly crib)**Product Code :**
BCIVP10011**Description :**

Baby Bassinet Deluxe (With poly crib) Frame work made of ERW tubular pipe of 2.5 cms. dia. Mounted on four 5 cms dia. castors. Removable CRIB of poly glass tiltable to trendelenburg and R-trend positions. Laminated bottom shelf with a spacious cabinet for storage. Provision for monitor stand and drip rod available. Pre-treated and epoxy powder coated

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,