



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Overbed Table (Adjustable With Geared Handle)

**Product Code :**  
NMIVP20031



**Description :**

Overbed Table (Adjustable With Geared Handle) Table base and frame work made of strong rectangular / square steel tubes mounted on four 5 cms dia. castors. Overall size of laminated top : 76L x 40W cms. adjustable by gear handle. Pre-treated and epoxy powder coated. Knockdown construction



**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,