



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Ambulance Foldaway Stretcher DX

Product Code :
NMSTR10003



Description :

Ambulance Foldaway Stretcher DX with wheels, outriggers and can be folded. It is made of high strength aluminum alloy and the water-tightness leather material. The back of this stretcher can be lifted, and the maximum angle is $60^{\circ} \pm 1^{\circ}$. It is characterized by its being fashionable structure, light-weighted, easy-carried, use-safely, reliable quality and easy for sterilization.

Dimensions

Product size (L x W x H): 185 x 50 x 24cm

Folding size (L x W x H): 92.5 x 50 x 11cm

Packing size (1pcs): 96 x 54 x 12cm Ambulance Foldaway Stretcher DX, Hospital Stretchers Exporters, Hospital Stretchers Manufacturers



Medical Hospital Equipments Manufacturers

Website: MedicalHospitalEquipmentsManufacturers.com