



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Overbed Table, Height Adjustable

**Product Code :**  
NMIVP20036



**Description :**

Overbed Table, Height Adjustable Overbed table with Laminated Wooden Top. Size: 750 ĩġ 1/2 350 ĩġ 1/2 820 / 1050 mm. Framework made of Stainless steel. Adjustable height from 820 - 1050 mm by gas spring. 4 Swivel castor



**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,