

**Product Name :**  
Comfy Hydraulic Stretcher Trolley, 5 Functions

**Product Code :**  
NMSTR10032



**Description :**

Comfy Hydraulic Stretcher Trolley, 5 Functions COMFY hydraulic stretcher trolley with adjustable back rest, knee rest, trendelenburg, reverse trendelenburg & adjustable height. Back Rest adjustable by gas spring, knee rest adjustable by super smooth crank, other functions adjustable by hydraulic mechanism. Size: 2082 x 860 x 608/923 mm. Framework made of mild steel epoxy coated finish with X-ray permeable Top. Collapsible side railings. Equipped with hydraulic pump. Equipped with I.V. Pole, drainage hook & bumper protection. 150mm dia noiseless castors with simultaneous braking system which locks / unlocks 2 castors with single pedal press. International quality structure and functionality to satisfy the needs of the best hospital and Operation theaters