

Product Name :
Hydraulic Stretcher Trolley, 2 Functions

Product Code :
NMSTR10034



Description :

Hydraulic Stretcher Trolley, 2 Functions Hydraulic stretcher trolley with adjustable Back rest & adjustable height. Back Rest adjustable by gas spring & Height adjustable by hydraulic mechanism. Size: 2100 $\frac{1}{2}$ 690 $\frac{1}{2}$ 460/850 mm. Framework made of mild steel epoxy coated finish with mesh mattress platform. Collapsible side railings. Equipped with bumper protection. 125mm dia noiseless castors with simultaneous braking system which locks / unlocks 2 castors with single pedal press. International quality structure and functionality to satisfy the needs of the best hospital and Operation theaters

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,