

**Product Name :**  
Stretcher Trolley**Product Code :**  
NMSTR10036**Description :**

Stretcher Trolley Size: 1950L x 650W x 680 / 900 H mm. Height adjustable from 600 - 950 H mm. Frame work made of stainless steel. S.S. Safety guard rails on each side. Mounted on four castor (all with individual brakes). Provision for IV rod. Supplied completely with S.S. side railings and IV rod

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers****Website:**Medical Hospital Equipments Manufacturers,