

Product Name :
Patient Transfer Stretcher, 2 Function

Product Code :
NMSTR10038



Description :

Patient Transfer Stretcher, 2 Function Size: 1940 x 680 x 610-920mm Height adjustable by crank mechanism (920 mm-610 mm) Parallel Mattress surface movement by crank to transfer patient to Bed or table. Horizontal movement range : 400mm Aluminum $\frac{1}{2}$ Alloy Guard rails, can be fixed up and down. 150mm Central $\frac{1}{2}$ locking castors, Equipped with separate directional wheel for easy maneuvering

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,