

Product Name :
Transfer Stretcher Cum Medical Trolley, Stainless Steel,
Adjustable Height

Product Code :
NMSTR10039



Description :

Transfer Stretcher Cum Medical Trolley, Stainless Steel, Adjustable Height Stainless Steel construction, for long life. Size: 3500 × 650 × 750 / 950 mm. Adjustable height from 750 - 950 mm by unobtrusive crank mechanism. Adjustable back section. The two stretcher base can move up / down independently and separately, and can be brought to the same level for smooth slide. Both base are locked securely before the bed slides. All wheels with individual brakes

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: MedicalHospitalEquipmentsManufacturers.com