



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Stretcher Trolley Aluminium With Adjustable Height,  
Adjustable Back rest

**Product Code :**  
NMSTR10045



**Description :**

Stretcher Trolley Aluminium With Adjustable Height, Adjustable Back rest Light Weight, easy clean, anodized aluminium alloy framework with adjustable back rest. Washable, water proof & anti-fungus fabric top. Size (LxWxH): 190x55x88cm, Size of Castors : 160mm. Height adjustable with crank mechanism. Supplied with I.V. Rod / Saline Rod. Range from bed surface to ground : 50cm. Self-Weight : 36 kg, Load Bearing: 159 kg



**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,