



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Electric Bed, Two Function For Home Care

Product Code :
NMHOPB00073



Description :

Electric Bed, Two Function For Home Care Two function electric bed for home care with electrically operated back rest: 0° to 70° $\frac{1}{2}$, knee rest: 0° to 40° $\frac{1}{2}$ $\frac{1}{2}$ $\frac{1}{2}$ $\frac{1}{2}$. Size: 2140mm(L)*1080mm(W)*500mm(H) All functions controlled with T-motion Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Equipped with battery backup. Epoxy coated mild steel frame work and 4 section perforated top. Wooden side panel Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy) Removable high quality beautiful wooden head panel and foot panel. 125mm dia noiseless castor, all with brakes Provision for I.V. Rod on all sides of the bed



Medical Hospital Equipments Manufacturers

Website: SITE