



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Sphygmomanometer Aneroid with Stethoscope

Product Code :
NMDIG10035



Description :

Sphygmomanometer Aneroid with Stethoscope Calibrated nylon cuff, stop pin manometer, 2-tube adult size bladder, Standard latex bulb, Special vinyl zipper case, Inflation Valve with spring, Standard end valve + Stethoscope



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,