



Email : sales@naugramedical.com

Phone: 9896600003

**Product Name :**  
Soiled Linen Trolley Conventional

**Product Code :**  
NMTOY10006



**Description :**

Conventional, Single Bag Tubular pipe frame work. Mounted on 75mm(3in) castor wheel. 480mm x 910mm x (19in x 36in) Bag. Canvas / Polyester fabric Bag.

Available in :

Epoxy powder coated.  
Stainless Steel.



**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)