



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Soiled Linen Trolley Conventional

Product Code :
NMTOY10006



Description :

Conventional, Single Bag Tubular pipe frame work. Mounted on 75mm(3in) castor wheel. 480mm x 910mm x (19in x 36in) Bag. Canvas / Polyester fabric Bag.

Available in :

Epoxy powder coated.
Stainless Steel.



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,