



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Professional Ophthalmoscope

**Product Code :**  
NMDIG10047



**Description :**

Professional Ophthalmoscope Halogen light for true tissue color and consistent, long-lasting illumination.  
Apertures: large and small spot, semicircle, and red-free filter, Correcting lenses range: +in 1D steps  
|1-4|6|8|10|15|20|0 & -in 1D steps |1-4|6|8|10|15|20 Bulb:2.5V 1.6W Halogen Lamp



**Medical Hospital Equipments Manufacturers**

Website: [SITE](http://SITE)