



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Fingertip Pulse Oximeter

Product Code :
NMDIG10049



Description :

Fingertip Pulse Oximeter Dependable Oxygen Level and Pulse Rate measurements Small & lightweight design Low power consumption Large LED Display of SpO2, PR, and Pulse Bar Auto shut-off with no signal for 5 seconds Suitable for adult, pediatric and neonate



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,