



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Fingertip Pulse Oximeter

Product Code :
NMDIG10050



Description :

Fingertip Pulse Oximeter Dependable Oxygen Level and Pulse Rate measurements LCD screen display SpO2, PR, and SpO2 waveform Small & lightweight design Low power consumption Auto shut-off with no signal for 5 seconds Build-in memory, with PC software analyzing Suitable for adult, pediatric and neonate



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,