

**Product Name :**  
Handheld Pulse Oximeter**Product Code :**  
NMDIG10051**Description :**

Handheld Pulse Oximeter Large LED Display of SpO2, PR, and Pulse Bar Durable, Compact and Lightweight Low power consumption Audible/Visual Alarms with Adjustable Alarm Limits Can Be Used for Spot-Check or Continuous Monitoring Suitable for Adult, Pediatric and Neonate

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**Website: [MedicalHospitalEquipmentsManufacturers.com](http://MedicalHospitalEquipmentsManufacturers.com)