

**Product Name :**  
Stretcher on Trolley**Product Code :**  
NMSTR10007**Description :**

Strong steel tubular frame work of trolley mounted on four 150 mm dia swiveling castors. Removable stretcher top. Made of aluminum sheet supported on tubular frame. Pre treated and Epoxy powder coated. Supplied with 35mm thick PU mattress.

Size: 2030 mm (L) x 560 mm (W) x 820 mm (H) (Approx.) Stretcher on Trolley, Hospital Stretchers Exporters, Hospital Stretchers Manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**Website: [MedicalHospitalEquipmentsManufacturers.com](http://MedicalHospitalEquipmentsManufacturers.com)