



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Ceiling Suspension Operation Theatre Light

Product Code :
NMOTL10002



Description :

28" DIA Satellite Dome with 7 Reflectors & 22" DIA Satellite Dome with 4 Reflectors, Fully spring balanced improved high intensity light for all major surgicals procedures, Both the domes are fitted with halogen bulbs of 24V -70W each, Provided with heat absorbing and colour correcting glass filters. Satellite dome fitted with a removable type sterilizable cum focussing handle for feather touch movements. Supplied with low voltage transformer

Light intensity :170000 Lux

Reflector : Diachronic

Colour Temp. : 4300K

Temp. Increased at head area : 2~4C

Finish : Oven Baked Ceiling Suspension Operation Theatre Light, Hospital OT Lights & Lamps Exporters, Hospital OT Lights & Lamps Manufacturers



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,