

Product Name :
Ultrasonic Nebulizers**Product Code :**
NMDIG10102**Description :**

Ultrasonic Nebulizers Durable, stretchable tube can produce a dense and tiny-particle mist. Power Voltage : AC220V \pm 10% Ultrasonic frequency : 1.7MHz \pm 10% Continuous working time : ? 4h Consumption power : ? 50W Max medication capacity for small cup : 150ml Max medication capacity for big cup : 350ml Max nebulizing rate: ? 3ml/min

NAUGRA®**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,