



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Surgical Operation Light

**Product Code :**  
NMOTL10003



**Description :**

28in DIA Satellite Dome with 7 Reflectors, Fully spring balanced improved high intensity light for all major surgicals procedures, Both the domes are fitted with halogen bulbs of 24V -70W each, Easy manoeuvrability by semi circular handle fitted outside the dome, Angular radial and axial movement easily manoeuvrable, Supplied with low voltage transformer

Light intensity :100000 Lux

Reflector : Diachronic

Colour Temp. : 4300K

Temp. Increased at head area : 2~4C

Finish : Oven Baked Surgical Operation Light, Hospital OT Lights & Lamps Exporters, Hospital OT Lights & Lamps Manufacturers



**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,