



Email : sales@naugramedical.com

Phone: 9896600003

**Product Name :**  
Digital Clinical Thermometer

**Product Code :**  
NMDIG10112



**Description :**

Digital Clinical Thermometer Numeric display - easy to read With Beep Alert Mercury Free (environmental friendly) Unbreakable High Precision Measurement Range: 32-42½°C & 89.6-109.4½°F Display Resolution: 0.1½°C & 0.1½°F Easy Replaceable Battery Auto power off 8 minutes Last measurement memory In plastic case Digital thermometer can be used orally, rectally & under armpit



**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,