



Email : [sales@naugramedical.com](mailto:sales@naugramedical.com)  
Phone: 9896600003

**Product Name :**  
Forehead Thermometer

**Product Code :**  
NMDIG10113



**Description :**

Forehead Thermometer Display range: 34.0~44.0oC/93.2~111.2oC Accuracy:  $\pm 0.3^{\circ}\text{C}/\pm 0.5^{\circ}\text{F}$  Min. scale: 0.1 Memory: 10 memories Measurement time (reference only): 1 second Beeper function Auto shut-off (1min after idle) Battery: 3V, Button battery(LR/CR 2032) Battery life: approx. 4,000 takes Dimensions: 140 x35x 35mm N.W.: 70g (including batteries) Safe and hygienic



**Medical Hospital Equipments Manufacturers**

Website: [SITE](http://SITE)