



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Digital EEG and Mapping System

Product Code :
NMDIG10116



Description :

Digital EEG and Mapping System Standard 24-channel(2400) amplifier The latest intelligent digital EEG amplifier, which integrates electrode box, amplifier, filter, A/D. Digital signal transmission enhances the anti jamming ability. 10/20 international standard system supports special electrodes, the combination of 13 kinds of standard leads can be selected (hardware) at real time or changed when uploading (software), support: A1, A2, A1+A2, A1A2. Resistance is displayed in the electrode box and tested in the acquisition course, so that the physician can adjust the electrode in time. USB interface, coupled with a laptop computer form a portable EEG system ECG and RESP Control flashlight automatically and manually The amplitude and frequency of EEG can be measured BEAM, power spectrum, and all kinds of trend graph The EEG database management, supporting MO or CD-RW archiving The special isolation transformer and optical cable ensure security Optional function Video monitor system Sleep module, ECG, EEG, SPO2, EMG, RESP Blue-tooth Communication



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,