



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Surgical Operation Light 20

Product Code :
NMOTL10004



Description :

Dome Diameter : 20in

Light intensity : 90000 Lux

Halogen Bulb : 24V - 150W

Reflector : Diachronic

Colour Temp. : 4300K

Temp. Increased at head area : 2~4C

Finish : Oven Baked Surgical Operation Light 20, Hospital OT Lights & Lamps Exporters, Hospital OT Lights & Lamps Manufacturers



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,