



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Mobile OT Light

Product Code :
NMOTL10005



Description :

Intensity: 70,000 Lux

Halogen Bulb: 24V - 70W

Reflector: Diachronic

Colour Temp.: 4300K

Temp. Increase at Head Area: 2-4C

Doma Dia: 22in

Finish: Oven Baked Mobile OT Light, Hospital OT Lights & Lamps Exporters, Hospital OT Lights & Lamps Manufacturers



Medical Hospital Equipments Manufacturers

Website: SITE,