

**Product Name :**  
Otoscope (Pin Connection System)**Product Code :**  
NMDIG10136**Description :**

Otoscope (Pin Connection System) 2.5V Standard illumination. Chrome plated handle, 4 x magnifying glasses. Reusable specula 2.5, 3.5, 4.5mm include. Connection for pneumatic Otoscope. Battery handle for 2x AA size batteries. Rheostat for adjustment of illumination. Supply : Black soft pouch or hard plastic case

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**Website: [MedicalHospitalEquipmentsManufacturers.com](http://MedicalHospitalEquipmentsManufacturers.com)