



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Professional Otoscope

Product Code :
NMDIG10138



Description :

Professional Otoscope Fiberoptic for cool, obstruction free light transmission Magnification:2.5x, Distal dia:2.5mm,4mm Bulb: 2.5V 1.6W halogen lamp Batteries: Two AA batteries or two "AA" (not included) NiMH rechargeable battery



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,