

Product Name :
Keratometer**Product Code :**
NMDIG10139**Description :**

Keratometer NET expertise in designing and manufacturing Keratometer that are used for measuring the curvature of the cornea for pre and post-surgical analysis and contact lens fitting. These keratometer are designed using latest technology and are thus reckoned for their ease of usage and accurate readings.

Technical Specifications Coincidence focusing system . Positive measurement to 0.02mm accuracy. 15x wide angle piece One position instrument. Measure both meridians without changing optical system. Simple vertical adjustable head to fit patient. Two way adjustable head and chin rests. Dual eye -level sighting system. Facilitates horizontal alignment. Positive fixation. Permit rapid measurement of central corneal area. Precession objectives for durability Achromatic Lenses. Measuring range: Radius of curvature, 9.4-6.4mm (in 0.05mm steps) and Diopter 36-52(in 0.025steps) Guarantee for one year. DIMENSIONS Height : 16" Weight : 7KGS Length : 19" Width : 7" Colour / Finish : Black Standard Accessories : Dust Cover, Two spares bulbs Optional Accessories at extra cost: Tripod stand/ Motorised instrument table.

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,